



Cold/warm compress thermPAD

TD-01-I.2.c-1

revision 1

date: 2023-05-24

amendment of: n/a

MANUFACTURER: ZARYS International Group sp. z o.o. sp.k.

ADDRESS: ul. Pod Borem 18, 41-808 Zabrze, Poland; T: +48 32 376 07 00-02,06 F: +48 32 370 38 08; www.zarys.com

CLASSIFICATION

Class I medical device compliant with Medical Devices Regulation (EU) 2017/745 (MDR)

INTENDED PURPOSE

- device used for the local cooling or heating of specific areas of the body to combat pain/neuralgia and reduce swelling and puffiness

FEATURES

- cold/warm gel compress made of non-toxic durable film with gel insert
- presence of a fabric cover to prevent external damage
- ability to store and return accumulated heat or cold for a long time
- remains flexible after freezing, allowing the compress to be easily shaped to the body area
- reusable
- non-sterile
- shelf life: 3 years



CATALOGUE NUMBER	DIMENSIONS	SALES PACKAGING		TRANSPORT PACKAGING	
		foil bag		carton	
KCZ7513	7,5 cm x 13 cm	1 pc.		100 x 1 pc.	
KCZ7552	7,5 cm x 52 cm	1 pc.		30 x 1 pc.	
KCZ1229	12 cm x 29 cm	1 pc.		25 x 1 pc.	
KCZ1314	13 cm x 14 cm	1 pc.		50 x 1 pc.	
KCZ1626	16 cm x 26 cm	1 pc.		20 x 1 pc.	
KCZ2138	21 cm x 38 cm	1 pc.		10 x 1 pc.	
KCZ3040	30 cm x 40 cm	1 pc.		8 x 1 pc.	