



Medline International France SAS  
5 rue Charles Lindbergh  
44110 Châteaubriant  
www.medline.com/en



## Technical sheet

# Flexible Enema Administration Set

### Reference(s)

**DYNDE70100**

### Description

Flexible Enema Administration Set contains:

- 1,5 l bag with tube, clamp and pre-lubricated tip
- 9ml liquid castile soap
- Fluid resistant under pad 42x45cm

### Main materials

Bag: Polyvinyl chloride (PVC)  
Tube: Polyvinyl chloride (PVC)  
Protective cap: Polyethylene resin (PE)  
Clamp: ABS  
Under pad: Cellulose / PE film

### Pack factor

48 units per carton

### Packaging

Carton box

### Manufacturer

Medline International France SAS  
EN ISO 13485

### Regulation

In compliance with Medical Devices Directive 93/42/EEC  
MDD Class: I



© 2017 Medline Industries Inc. or one of its subsidiaries.  
All rights reserved.

LIT-Enema Kit  
DYNDE70100-  
EN01

01.01.INT.12