

4560-164 Guilhufe, Penafiel - Portugal

Ref: 4881-502

Product Specification

Brand:



DRQ 14/21 REV.00 2020/09/08

Description: Utility drape

Final dimensions:	45 cm x75 cm
Composition:	Polyethylene film,SB polypropylene
	nonwoven.
Size variation:	±10 %
Latex:	Free
PVC:	Free
Re use:	No
Colour:	green
Sterilized by:	Ethylene Oxide
Shelf life:	60 months
CFU/g	N.A.
Class.EN 13795-1:2019	STANDARD PERFORMANCE
Resist. Micro Penetr. Wet	>2,8 IB
Cleanliness- Microbial	N.A. CFU/100cm^2
Clean particulate matter	<3,5 IPM
Linting	<4,0 log10lint count
Resist. Liq.Penetration	>30 cmH2O
Dry Bursting Strength	>40 kPa
Wet bursting strength	>40 kPa
Dry Tensile Strength	MD>40/CD>20 N/5cm
Wet Tensile Strength	MD>40/CD>20 N/5 cm





Nr and Revision Date:	16 2020/12/18		
CE Class:	I Steril		
CE mark:	CE1639		
Compliance with:	MDD 93/42/EEC amended by 2007/47/EC		
	EN_13795_1_2019 ; EN_ISO_13485.2016 ;	EN_ISO_13485.2016_AC18 ; EN 556-1.2001_AC.2006 ;	
	ISO_11135.2014; ISO_10993-1.2018; EN_	_ISO_11607_1 2017; EN_ISO_11607-2.2017;	
	ISO_TS_16775.2014; ASTM_F1980_16; ISO_14644-1:2015; ISO_14698_1_2003;		
Annlicable norms:	FN ISO 14971.2012 : FN ISO 15223-1.201	16 : FN 1041.2008 A1.2013 :	

Level		Туре	Dimension	Stamp	Labels
1	1	Flexible peel-open pack	134x200mm		
2	90	Folding box	216x151x308mn	n	labeled
3	900	Carton box	780x225x635mn	n	labeled

	The device is intended to be used on a single patient during one single procedure to prevent risks from	
Important Notes:	contamination, including possible infections, and/or functional damage.	

	Store in a dry place at room temperature and protect from sunlight. The product has a good stability and therefore, normal conditions of medical devices storage shall apply	
	No special requirements. The waste must be disposed according to local legal regulations for the disposable of medical waste.	
Intended use	As covering drape to prevent transfer of infective agents in surgical context.	

Environment - Packaging waste	Weight (gr) of packing material used in 100 pcs	634
-------------------------------	---	-----

ISSUED BY: Quality Management Department, 2020/12/18

CONFIDENTIAL: The information contained herein is confidential and belongs to Bastos Viegas, s.a. It is provided for specific purposes and in connection with established agreements. It is not permitted disclosure, reproduction, copying or distribution, in whole or in part, for other purposes or to any third parties.